

## Alcian Blue (pH 2.5) Stain Kit

**Description:** The Alcian Blue (pH 2.5) Stain Kit is intended for use in the histological visualization of sulfated and carboxylated acid mucopolysaccharides and sulfated and carboxylated sialomucins (glycoproteins).

Acidic Sulfated Mucosubstances: Blue  
Hyaluronic Acid: Blue  
Sialomucins: Blue  
Nuclei: Red  
Background: Pink

**Uses/Limitations:** For In-Vitro Diagnostic use only.  
Histological applications.  
Do not use past expiration date.  
Use caution when handling these reagents.

**Control Tissue:** Small Intestine  
Appendix  
Colon


### Availability/Contents:

<u>Kit Contents</u>	<u>Volume</u>	<u>Storage</u>
Alcian Blue, pH 2.5	250 ml	15-30°C
Nuclear Fast Red Solution	250 ml	15-30°C
Acetic Acid Solution	500 ml	15-30°C


**Precautions:** This product is a single-use, non-sterile, in vitro diagnostic device.  
Avoid contact with skin and eyes.  
Harmful if swallowed.  
Follow all Federal, State, and local regulations regarding disposal.  
Use in chemical fume hood whenever possible.

### Procedure (Standard):

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Incubate slide in Acetic Acid solution for 3 minutes.
3. Incubate slide in Alcian Blue (pH 2.5) solution for 30 minutes at room temperature or 15 minutes at 37° C.
4. If desired, rinse slide briefly in Acetic Acid solution to remove excess Alcian Blue.
5. Rinse for 2 minutes in running tap water followed by 2 changes of distilled water.
6. Stain slide in Nuclear Fast Red Solution for 5 minutes.
7. Rinse for 2 minutes in running tap water followed by 2 changes of distilled water.

Storage: 15° C  30° C

**Store All Components At  
Room Temperature.**



Diagnostic BioSystems, Inc.  
1020 Serpentine Lane, Suite 114  
Pleasanton, CA 94566  
U.S.A.

**C V  
P**

EmergoEurope  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

# Instructions For Use KT 003-IFU

Rev. Date: 2/17/2022

**Revision: 5**

Page 2 of 2

6616 Owens Drive, Pleasanton, CA 94588 U.S.A. - Tel. (925) 484-3350 - Fax (925) 484-3390 - [www.dbiosys.com](http://www.dbiosys.com)

8. Dehydrate through graded alcohols.
9. Clear, and mount in synthetic resin.

## References:

1. Lillie, R.D. 1977, H.J. Conn's Biological Stains, 9<sup>th</sup> Edition. Williams & Wilkins, Baltimore. Pages 452-455.
2. Sheenan, D.C., Hrapchak, B.B. Theory and Practice of Histotechnology, 2<sup>nd</sup> Edition. Battelle Press, Columbus, OH. Pages 172-173.
3. Churukian, C.J., 1989, Manual of Special Stains Laboratory, 4<sup>th</sup> Edition. University of Rochester, Rochester, New York. Pages 55-56.
4. Carson, F.L., 1996, Histotechnology; A Self-Instructional Text, 2<sup>nd</sup> Edition. ASCP Press, Chicago, IL. Pages 117-121.

Storage: 15° C



30° C

**Store All Components At  
Room Temperature.**



Diagnostic BioSystems, Inc.  
1020 Serpentine Lane, Suite 114  
Pleasanton, CA 94566  
U.S.A.

**C V  
P**

EmergoEurope  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands