

Fite's Stain Kit

(For Leprosy and Nocardia)

Description: The Fite's Stain Kit (For Leprosy and Nocardia) is intended for use in the histological visualization of mycobacterium Lepra bacillus and Nocardia. This kit may be used on formalin-fixed, paraffin-embedded sections.

Lepra bacillus:	Red
Nocardia:	Red
Background:	Blue

Uses/Limitations: For In-Vitro Diagnostic use only.
Histological applications.
Do not use past expiration date.
Use caution when handling these reagents.

Control Tissue: Any well fixed paraffin embedded tissue.


Availability/Contents:

<u>Kit Contents</u>	<u>Volume</u>	<u>Storage</u>
Xylene-Peanut Oil Solution	125 ml	15-30°C
Carbol Fuchsin Solution	125 ml	15-30°C
Acid Alcohol Solution (1%)	500 ml	15-30°C
Methylene Blue Solution	125 ml	15-30°C


Precautions: This product is a single-use, non-sterile, in vitro diagnostic device.
Keep away from open flame.
Avoid contact with skin and eyes.
Harmful if swallowed.
Follow all Federal, State, and local regulations regarding disposal.
Use in chemical fume hood whenever possible.

Leptra bacillus Procedure (Standard):

1. Deparaffinize sections in 2 changes of Xylene-Peanut Oil Solution for 12 minutes each.
2. Air dry slide for 15 minutes. Do not remove oil film. Remaining film prevents de-staining of Lepra bacillus during differentiation.
3. Rinse slide in several changes of distilled water.
4. Incubate slide in Carbol Fuchsin Solution for 15 minutes.
5. Rinse slide in several changes of distilled water.
6. Differentiate section in Acid Alcohol Solution (1%) until background is pink.
7. Rinse slide in distilled water and check by microscope for correct differentiation.

Storage: 15° C  30° C

**Store All Components At
Room Temperature.**



Diagnostic BioSystems, Inc.
6616 Owens Drive
Pleasanton, CA 94588
U.S.A.

C V
P
Emergo Europe
Prinsessegracht 20
2514 AP The Hague, The Netherlands

6616 Owens Drive, Pleasanton, CA 94588 U.S.A. - Tel. (925) 484-3350 - Fax (925) 484-3390 - www.dbiosys.com

8. Rinse in running tap water for 1 minute followed by 1 rinse in distilled water.
9. Dip slide 2-3 times in Methylene Blue Solution.
10. Dip slide quickly in distilled water and check by microscope for correct staining.
11. Air dry slide at room temperature.
12. Dip slide several times in Xylene or Xylene Substitute.
13. Mount in synthetic resin.

Nocardia Procedure:

Preparation of Reagents Prior to Beginning:


1. Prepare **Diluted Acid Alcohol Solution** by mixing 25ml of Acid Alcohol Solution (1%) with 25ml of Distilled Water.

Procedure:


1. Deparaffinize sections in 2 changes of Xylene-Peanut Oil Solution for 12 minutes each.
2. Air dry slide for 15 minutes. Do not remove oil film. Remaining film prevents de-staining of Lepra bacillus during differentiation.
3. Rinse slide in several changes of distilled water.
4. Incubate slide in Carbol Fuchsin Solution for 15 minutes.
5. Rinse slide in several changes of distilled water.
6. Dip slide 2-3 times in Diluted Acid Alcohol Solution.
7. Rinse slide in distilled water and check by microscope for correct differentiation. Avoid decolorizing the Nocardia organism.
8. Rinse in running tap water for 1 minute followed by 1 rinse in distilled water.
9. Dip slide 2-3 times in Methylene Blue Solution.
10. Dip slide quickly in distilled water and check by microscope for correct staining.
11. Air dry slide at room temperature.
12. Dip slide several times in Xylene or Xylene Substitute.
13. Mount in synthetic resin.

References:

1. Fite, G.L. Cambre, P.J., and Turner, M.H. Procedure for demonstrating lepra bacilli I paraffin sections. Arch. Pathol., Vol. 43, pages 624-25, 1947.
2. Mallory, Pathological Technique; page 275.
3. Clarke, G., et al. Staining Procedures, 4th Edition, Williams& Wilkins, page 447, 1981.

Storage: 15° C  30° C

**Store All Components At
Room Temperature.**

 Diagnostic BioSystems, Inc.
6616 Owens Drive
Pleasanton, CA 94588
U.S.A.

C V
P
Emergo Europe
Prinsessegracht 20
2514 AP The Hague, The Netherlands

Instructions For Use KT013-IFU


Rev. Date: Dec. 12, 2018

Revision: 6

Page 3 of 3

6616 Owens Drive, Pleasanton, CA 94588 U.S.A. - Tel. (925) 484-3350 - Fax (925) 484-3390 - www.dbiosys.com

4. Crowder, C., Taylor, HW., Modified Fite Stain for Demonstration of Mycobacterium Species in Tissue Sections; Journal of Histotechnology; Volume 19; 2: pages 133-134.

Storage: 15° C  30° C

**Store All Components At
Room Temperature.**



Diagnostic BioSystems, Inc.
6616 Owens Drive
Pleasanton, CA 94588
U.S.A.

C V
P

Emergo Europe
Prinsessegracht 20
2514 AP The Hague, The Netherlands