



Instructions For Use

KT 008-IFU

Document #: DS-3003-D

Release Date: 07/07/2025

Amyloid Stain Kit (Congo Red)

Description and Principle

The Amyloid Stain Kit (Congo Red) is intended for use in the histological visualization of amyloid in tissue sections. Examination under a polarizing microscope results in green birefringence of amyloid. Congo Red binds amyloid by hydrogen bonds with amyloid fibers in a highly oriented linear manner. Linear bound amyloid will exhibit green birefringence under polarized light.

Expected Results

Amyloid:	Red to Pink
Erythrocytes:	Light Orange
Eosinophil Granules:	Orange to Red
Nuclei:	Blue

Kit Contents

	Volume	Storage
1. Congo Red Solution	500 ml	15-30°C
2. Hematoxylin	500 ml	15-30°C
3. Bluing Reagent	500 ml	15-30°C

Suggested Controls (not provided)

Freshly cut FFPE sections containing amyloid.
It has been reported that birefringence of pre-cut sections weakens over time.
Cut sections 8-10 microns for optimal visualization of amyloid deposit birefringence.

Uses/Limitations

For In-Vitro Diagnostic use only.
Do not use if reagents become cloudy or precipitate.
Do not use past expiration date.
Use caution when handling reagents.
Non-Sterile
Intended for FFPE sections cut at 5-10µm.
This procedure has not been optimized for frozen sections.
Frozen sections may require protocol modification.

Storage

Store kit and all components at room temperature (15-30°C).

Safety and Precautions


Please see current Safety Data Sheets (SDS) for this product and components GHS classification, pictograms, and full hazard/precautionary statements. If there is any serious incident that has occurred in relation to the device, please contact the manufacturer: Diagnostic BioSystems Technical Support at (925) 484-3350, extension 2 or techsupport@dbiosys.com. If required, please report to the Competent Authority of the Member State in which the user and/or patient is established.

Procedure

1. Deparaffinize sections and hydrate to distilled water.
2. Stain slide with Hematoxylin for 5 minutes.
3. Rinse slide in tap water.
4. Incubate slide in Bluing Reagent for 30 seconds.

5. Rinse slide in distilled water.

6. Briefly rinse slide in absolute alcohol for 3-5 seconds.



Amyloidosis in Canine Kidney demonstrated with Congo Red and viewed under A) Brightfield, and B) Polarized Light

7. Pour Congo Red Solution into a staining jar and place slide(s) in for 20 minutes. Ensure there is enough stain to completely cover tissue. Staining in smaller volumes of Congo Red Solution by pipetting or pouring on horizontal slides may cause non-specific staining.

8. Rinse excess stain off slide with absolute alcohol.

9. Dehydrate slide in absolute alcohol for another 30 seconds.

10. Clear and mount in synthetic resin.



Diagnostic BioSystems
6616 Owens Drive
Pleasanton, CA, 94588
Tel: (925) 484 3350
www.dbiosys.com



CH REP
MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland



MedEnvoy Global B.V.
Prinses Margrietplantsoen 33 - Suite 123
2595 AM The Hague
The Netherlands

References


1. Puchtler, H, et al: On the binding of Congo Red amyloid. J. Histochem. Cytochem. Vol. 10: pages 355-363, 1962.
2. Eastwood, H. & Cole, K.R., Staining of amyloid in buffered Congo Red in 50% ethanol. Stain Technology. Vol. 46: pages 208-209, 1971.
3. Carson, F.L., Histotechnology: A Self-Instructional Text, 2nd Edition. ASCP Press, Chicago, IL. Pages 117-121, 1996.
4. Churukian, C., Improved Puchtler's Congo Red method. J. of Histotechnology. Vol. 23: pages 139-141, 2000.



Diagnostic BioSystems
6616 Owens Drive
Pleasanton, CA, 94588
Tel: (925) 484 3350
www.dbiosys.com



CH	REP	
MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland		

EC	REP		MedEnvoy Global B.V. Prinses Margrietplantsoen 33 - Suite 123 2595 AM The Hague The Netherlands
----	-----	---	--