

# SITVue™

## Three Step DAB Detection System

**SPEED  
MATTERS**  
RESULT IN




### Anti Mouse/Rabbit SITVue™ HRP/DAB Detection System


SITVue™/DAB Detection System is a powerful Intensification system that significantly enhances chromogenic signals. This system can be integrated into standard immunohistochemistry (IHC) staining methods.

SITVue™/DAB Detection System is an enzyme-mediated reaction that utilizes horseradish peroxidase (HRP) to catalyze the deposition of two separate Linkers applied sequentially onto tissue sections or cell preparation. The deposited Linkers can be detected with a streptavidin-peroxidase conjugate followed by a reaction with a peroxidase substrate/chromogen solution such as diaminobenzidine (DAB).


The SITVue™/DAB Detection system results in a significant increase in sensitivity compared to standard IHC detection methods, while maintaining similar specificity.



SITVue™ /DAB Detection System




Ki67




Cytokeratin Pan

Two Steps Detection System



Ki67



Cytokeratin Pan

### Ordering Information

#	Cat. No	Volume
1	SIT-1000D	1000 Tests
2	SIT-100D	100 Tests




**Diagnostic  
BioSystems**

Step	Staining Procedure	Two Steps Detection system Incubation Time	SITVue™ /DAB Detection System Incubation Time
Peroxidase Block	A. Incubate slides in Peroxidase Block.	5 min.	1min.
	B. Rinse slides with Immuno Wash Buffer three (3) times, for 1 min. each time.	3 x 1 min.	
Pre-Blocking (optional)	A. Add 2 drops (100 µL) or enough volume of Pre-Blocking Solution to cover the tissue section.	10 min.	N/A
	B. Drain or blot off solution. Do not rinse.		
Primary Antibody	A. Incubate with Primary Antibody, prepared according to the manufacturer's recommended protocol at the desired concentration. Concentrated Primary Antibodies may be diluted using Primary Antibody Diluent.	30 – 60 min.	5 min.
	B. Wash slides with 3 changes of Immuno Wash Buffer.	3 x 1 min.	3 x 10 second.
Linker 1	A. Apply the Linker 1 and incubate.	10 min.	5 min.
	B. Rinse slides with 3 changes of Immuno Wash Buffer.	3 x 1 min.	3 x 10 second.
Linker 2	A. Incubate the tissue with Linker 2.	N/A	5 min.
	B. Wash slides with 3 changes of Immuno Wash Buffer.		3 x 10 second.
Tracer	A. Incubate the tissue with Tracer.	10 min.	5 min.
	B. Wash slides with 3 changes of Immuno Wash Buffer.	3 x 1 min.	3 x 10 second.
Stable DAB/Plus	A. Prepare the Stable DAB/Plus substrate working solution (see above).		
	B. Incubate tissue with prepared Stable DAB/Plus substrate solution. Monitor level of staining to determine optimal time of incubation.	5 – 10 min.	5 min.
	C. Rinse slides with 3 changes of water.	3 x 1 min.	3 x 10 second.
Counterstain	A. Incubate tissue with Counterstain (Hematoxylin K097), according to manufacturer's recommendation or standard laboratory protocol.	~1 min.	~1 min.
	B. Wash slides with water 3 times, followed by 1 time in Immuno Wash Buffer, then 1 time in water.	3 x 1 min. H2O 1 x 1 min. Buffer 1 x 1 min. H2O	3 x 10 second. H2O 1 x 10 second. Buffer 1 x 10second. H2O
<b>Total staining time</b>		<b>90-125 min</b>	<b>&lt; 30 min</b>
Dehydrate & Coverslip	A. Dehydrate tissues through graded ethanol series, followed by xylene series. B. Apply coverslips with permanent mounting medium.		



#### United States (Headquarters)

6616 Owens Drive, Pleasanton, CA 94588, USA

Toll-Free (USA)  : +1 (844)856-3350, customersupport@dbiosys.com



Contact your distributor  
for more information  
and visit us at [www.dbiosys.com](http://www.dbiosys.com)

#### Europe

Germany  
Lake Constance

Cell: +49 7557 929 3915  
[wolfgang.vogel@dbiosys.com](mailto:wolfgang.vogel@dbiosys.com)

#### Mexico

Avenida Aztecas # 95, Casa 4  
Colonia Pueblo de los Reyes  
Alcaldia Coyoacan  
CDMX, C.P. 04330

Cell: +525535709693  
[a.herrera@dbiosys.com](mailto:a.herrera@dbiosys.com)

#### India

612, Eden Square  
St. John's Road  
Secunderabad – 500003

Cell: +91 9958293222  
[anandam@dbiosys.com](mailto:anandam@dbiosys.com)

#### Canada

32 Lemsford Drive  
Markham, ON L3S 4H4

Cell: +1 416 219 2035  
[hosna.mujadidi@dbiosys.com](mailto:hosna.mujadidi@dbiosys.com)