

Concentrated Alkaline Phosphatase-conjugated Streptavidin

Catalog No.:	C 006
Format:	1mL cherry red solution containing Alkaline Phosphatase-conjugated Steptavidin.
Intended Use:	This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy. Clinical interpretation of staining results should be accompanied by histological studies with proper controls. Patients' clinical histories and other relevant diagnostic tests should be utilized by a qualified person(s) when evaluating and interpreting results.
Suggested Dilution:	Dilute 1:100 to 1:200 for immunohistochemistry.
Buffer:	0.01M Sodium Phosphate Buffer, pH 8.0.
Stabilizer:	5 mg/mL Bovine Serum Albumin
Preservative:	0.1% Sodium azide
Precautions:	Not intended for household use. Do not take internally. May be harmful or fatal if swallowed. Avoid contact with eyes, skin and clothing.

IVD: For In Vitro Diagnostic Use

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

DBS

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: www.dbiosys.com e-mail: customersupport@dbiosys.com