

## Goat Polyclonal Antibody to Hepatitis B Surface Antigen (HBsAg)

<b>Catalog No.:</b>	GP 001, GP 001-05
<b>Applications:</b>	Analyte Specific Reagent. Analytical and Performance characteristics are not established. Optimal working conditions must be determined by the end user.
<b>Immunogen:</b>	Purified surface antigen of subtype ad/ay.
<b>Host:</b>	Goat
<b>Format:</b>	Partially purified goat antiserum containing sodium azide as a preservative.
<b>Specificity:</b>	This antigen reacts with liver cells of patients infected with Hepatitis type B virus and with hepatocellular carcinoma and cirrhosis.
<b>Positive Control:</b>	Infected liver
<b>Cellular Localization:</b>	Cytoplasmic
<b>Storage:</b>	Store at 2-8°C. Do not use beyond the expiration date stated on the label.

**ASR: Analyte Specific Reagent**

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

**DBS**

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: [www.dbiosys.com](http://www.dbiosys.com) e-mail: [customersupport@dbiosys.com](mailto:customersupport@dbiosys.com)