

# ALP Label

<b>Catalog No.:</b>	M 009
<b>Product:</b>	Streptavidin-conjugated Alkaline Phosphatase
<b>Applications:</b>	For immunohistochemical staining of tissue sections. Ready to Use. Incubate for 20-30 minutes at room temperature.
<b>Format:</b>	10mL clear cherry red colored solution of streptavidin-conjugated alkaline phosphatase.
<b>Buffer:</b>	0.01M Tris buffer saline Buffer, pH 7.6.
<b>Stabilizer:</b>	1mg/mL Bovine Serum Albumin
<b>Preservative:</b>	0.05% Proclin 300
<b>Storage:</b>	Store at 2-8°C. Do not use beyond the expiration date stated on the label.
<b>Precautions:</b>	<ul style="list-style-type: none"><li>(i) Do not take internally. May be harmful or fatal if swallowed. Avoid contact with eyes, skin, and clothing.</li><li>(ii) This reagent contains Sodium Azide. Follow instructions provided by local authorities for disposal. If disposed in the sink, flush the drain pipe with water to avoid a reaction of Sodium Azide with the plumbing system.</li><li>(iii) Interpretation of the results will be the sole responsibility of the user.</li></ul>

## IVD: For In Vitro Diagnostic Use

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

**DBS**

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: [www.dbiosys.com](http://www.dbiosys.com) e-mail: [customersupport@dbiosys.com](mailto:customersupport@dbiosys.com)