

Goat Polyclonal Antibody to Hepatitis B Surface Antigen (HBsAg)

Catalog No.:	PDR 031
Applications:	Analyte Specific Reagent. Analytical and Performance characteristics are not established. Optimal working conditions must be determined by the end user.
Immunogen:	Purified surface antigen of subtype ad/ay.
Host:	Goat
Specificity:	This antigen reacts with liver cells of patients infected with Hepatitis type B virus, hepatocellular carcinoma, and cirrhosis.
Positive Control:	Infected liver
Cellular Localization:	Cytoplasmic
Storage:	Store at 2-8°C. Do not use beyond the expiration date stated on the label.

ASR: Analyte Specific Reagent

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

DBS

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: www.dbiosys.com e-mail: customersupport@dbiosys.com