

Immunohistochemistry Humidity Chamber


Diagnostic BioSystems Slide Master365™ is a compact and efficient platform that facilitates Humidity required for IHC staining. Appropriate maintenance of humidity along with flat platform is of extreme importance for consistent immunohistochemistry staining.

Capabilities

- ▶ Ensures flat platform.
- ▶ Maintains humidity.
- ▶ Ease to rinse and avoids majority of individual slides handling during staining.
- ▶ 40 slide staining capacity.



Specifications	40 Slides
Color	Black and White
Weight	2.5 kg
Dimensions	470 x 370 x 50 mm
Slide Capacity	40 to 50 Slides
Optimal Run Temperature	Room Temperature



United States (Headquarter)

6616 Owens Drive, Pleasanton, CA 94588, USA
Toll-Free (USA) 🇺🇸 : +1 (844) 856-3350
customersupport@dbiosys.com



www.dbiosys.com