



Explore...
Diagnostic BioSystems
Ancillary Reagents

Diagnostic BioSystems offers a wide range of ancillary reagents for varied immunohistochemistry applications. Our ancillary products include antigen retrieval solutions, buffers, mounting media, and other accessory products, designed to provide solutions to your IHC challenges. All of our products are manufactured in accordance with ISO 9001:2008 and FDA cGMP guidelines.

Affordable Solutions For Your IHC Staining Challenges

Diagnostic BioSystems Ancillary Reagents



Pens



Buffers



Antigen Retrieval Solutions



Mounting Media

Images not to scale

Antigen Retrieval Solutions

Diagnostic BioSystems offers a wide range of products for both heat induced and enzyme based antigen retrieval.

Product	Volume	Catalog No:
HEAT INDUCED ANTIGEN RETRIEVAL		
10X Antigen Retrieval Buffer, pH 10.0	500ml	K 044
10X Citrate Buffer, pH 6.0	500ml	K 035
10X Citrate Buffer, pH 7.0	500ml	K 036
10X Citrate Buffer, pH 8.0	500ml	K 037
10X EDTA Buffer, pH 8.0	500ml	K 038
10X Tris-EDTA Buffer, pH 9.0	500ml	K 043
ENZYME PRETREATMENT		
25X Proteinase K	2ml	K 030
Antigen Unmasker System	20 tablets	DSI5030
HistoZyme	15ml, 50ml	K 046, K 046-50
Solution Pepsin	50ml	K 009

Diagnostic BioSystems is dedicated to providing innovative solutions to your immunohistochemistry needs. For more information about our products and solutions, visit us at: www.dbiosys.com OR on Facebook at: www.facebook.com/DiagnosticBioSystems.

Mounting Media

We offer a broad selection of easy to use permanent and aqueous slide mounting media for tissue and cell visualization.

Product	Volume	Catalog No:
CC/Mount™	30ml, 200ml	K 002, K 002-200
Flouromount™	25ml , 200ml	K 024, K 024-200
Flouromount/Plus™	25ml, 200ml	K 048, K 048-200

Buffers

Our wash buffer is specifically designed to remove loosely bound antibodies effectively and efficiently to provide a cleaner background.

Product	Volume	Catalog No:
10X Immuno Wash Buffer	1000ml	K 005
Phosphate Buffered Saline Tablets	100 tablets	DMR E404-01

Blocking Agents

We offer blocking agents which eliminate non-specific binding and background noise

Product	Volume	Catalog No:
Background Blocker™	100ml	K 023
Peroxidase Block	50ml	K 033
Tissue Primer	100ml	K 054

Pens

Our PAP Pen's prevent the waste of valuable reagents by creating a water repellent barrier.

Product	Quantity	Catalog No:
Elite Mini PAP Pen	1 pen	K 042
Elite PAP Pen	1 pen	K 039

Accessory Products

Product	Volume	Catalog No:
10X Peroxidase Diluent I	100ml	K 019
10X Peroxidase Diluent II	100ml	K 020
10X Secondary Antibody Diluent	100ml	K 022
Alkaline Phosphatase Diluent	1000ml	K 028
Alkaline Phosphatase Enhancer	250ml	K 011
Broad Spectrum Negative Control	6ml	K 053-06
	25ml	K 053-25
	100ml	K 053-100
Peroxidase Block	50ml	K 033
Primary Antibody Diluent	250ml	K 004
TissueBond™ Reagent	7ml	K 013



**Diagnostic
BioSystems**

Our Innovation. Your Solution.



**Diagnostic
BioSystems**

United States (Headquarters)

6616 Owens Drive, Pleasanton, CA 94588, USA

Toll Free: (888) 896-3350, customersupport@dbiosys.com

Europe

Germany
Lake Constance

Cell: +49 7557 929 3915
wolfgang.vogel@dbiosys.com

Mexico

Avenida Aztecas # 95, Casa 4
Colonia Pueblo de los Reyes
Alcaldia Coyoacan
CDMX, C.P. 04330

Cell: +525535709693
a.herrera@dbiosys.com

India

612, Eden Square
St. John's Road
Secunderabad – 500003

Cell: +91 9958293222
anandam@dbiosys.com

Canada

32 Lemsford Drive
Markham, ON L3S 4H4

Cell: +1 416 219 2035
hosna.mujadidi@dbiosys.com



Contact your distributor
for more information
and visit us at www.dbiosys.com