

## Immunohistochemistry Humidity Chamber

Diagnostic BioSystems Slide Master365™ is a compact and efficient platform that facilitates Humidity required for IHC staining. Appropriate maintenance of humidity along with flat platform is of extreme importance for consistent immunohistochemistry staining.

### Capabilities

- ▶ Ensures flat platform.
- ▶ Maintains humidity.
- ▶ Ease to rinse and avoids majority of individual slides handling during staining.
- ▶ 40 slide staining capacity.



The Innovative tissue staining humidity chamber



Slotted drain holes work to maintain humidity



Rinse in one easy motion

Specifications	40 Slides
Color	Black and White
Weight	2.5 kg
Dimensions	470 x 370 x 50 mm
Slide Capacity	40 to 50 Slides
Optimal Run Temperature	Room Temperature

#### Headquarters

6616 Owens Drive  
Pleasanton, CA 94588, USA  
Toll Free: (888) 896-3350  
customersupport@dbiosys.com



[www.dbiosys.com](http://www.dbiosys.com)