

A20 (A-12) Antibody


Mouse Monoclonal Antibody

A20, also referred to as TNF- α -induced protein 3 (TNFAIP3), is cytokine-inducible protein that functions to inhibit apoptosis and activate NF- κ B. It was first identified as a TNF- α inducible primary response gene in human umbilical vein endothelial cells, and encodes a 790-amino acid protein containing seven Cys2/Cys2-zinc finger motifs. Expression of A20 can confer resistance to apoptosis and NF- κ B activation triggered by these signals, probably through interference with TRAF (TNF receptor associated factor) family members, and interaction with the NF- κ B inhibiting protein ABIN.


Antibody Name	Cat. No	Clone	Type	Localization	Format	Regulatory
A20	Mob598/ PDM598	A-12	Mouse Monoclonal	Cytoplasmic	Concentrate & Prediluted	IVD

Why A20?

- Constitutive expression of A20 is observed in lymphoid tissues
- Haploinsufficiency of TNFAIP3 (A20) by a germline heterozygous mutation leads to the ALPS phenotype.
- A20 are frequent events in Diffuse large B-cell lymphoma (DLBCL)



Formalin fixed paraffin
embedded human Lung Ca.
stained with A20 25x



Formalin fixed paraffin
embedded human Lung Ca.
stained with A20 40x

Diagnostic BioSystems

6616 Owens Drive
Pleasanton, CA 94588, USA
T: 925-484-3350
www.dbiosys.com

Toll-Free (USA)  : +1 (844) 856-3350
customersupport@dbiosys.com

Reference:

1. Beyaert, R. et al. (2000) Biochem Pharmacol 60, 1143-51.
2. Song, H.Y. et al. (1996) Proc Natl Acad Sci USA 93, 6721-5.
3. Heyninck, et al. (1999) FEBS Lett 442, 147-50.