



CD56 Antibody

Mouse Monoclonal Antibody

CD56 is extremely useful in the diagnosis of papillary thyroid carcinoma (PTC), including the follicular variant, and to distinguish it from other follicular cell-derived thyroid tumors/lesions. Used to differentiate plasma cells in myeloma (CD56+) from reactive plasmacytosis or MGUS (CD56-), also to predicts poor prognosis for cutaneous lymphoproliferative disorders other than cutaneous T cell lymphoma.

Antibody Name	Cat. No	Clone	Type	Localization	Format	Regulatory
CD56	Mob261/ PDM110	123C3.D5	Mouse Monoclonal	Cell membrane	Concentrated & Prediluted	IVD

Why CD56




- CD56 is the prototypic marker of NK cells (80-90%), and NK/T-cell lymphomas
- To differential diagnosis of Myeloid leukemic cells (CD56+) with normal myeloid cells (CD56-)
- To differential non-papillary thyroid tumors (CD56+) with papillary thyroid carcinoma (CD56-)
- To differential Schwannoma (CD56+) with neurofibroma (CD56-)




Nordic Immunohistochemical Quality Control Assessment of CD56

Assessment 61- individual results: **Good**



Fomalin fixed paraffin
embedded human Colon Ca.
stained with CD56 10x



Fomalin fixed paraffin
embedded human Colon Ca.
stained with CD56 40x

Diagnostic BioSystems

6616 Owens Drive
Pleasanton, CA 94588, USA
T: 925-484-3350
www.dbiosys.com

Toll-Free (USA)  : +1 (844) 856-3350
customersupport@dbiosys.com