

Ensures high quality IHC staining

The Diagnostic BioSystems



Antigen Retrieval System
For Crisp and Clean staining

Features

- Up to 5 cycles per run with individual time and temperature settings
- Real time temperature and time monitoring
- Low reagent consumption
- Minimal reagent evaporation
- Uniform retrieval, even in over fixed tissues
- Multiple retrieval buffers in a single run
- IHC, ISH, and FISH suitability
- Add on slides in-between the cycles
- Temperature ranges from ambient to 115 °C
- High throughput (96 slides in less than 30 minutes)



Cat.No: AR365



Visit us at www.dbiosys.com for more details.

*Country specific power cord need to be procured locally.

Montage Opus365™

Antigen Retrieval System



Each **Montage Opus365™** unit includes the following



Antigen Retrieval Tanks (4)



24 Slide Rack (4)



Metal Stand (1)



Montage 10X EDTA pH 8.0 AR Solution (1)



Montage 10X Citrate pH 6.0 AR Solution (1)



**Diagnostic
BioSystems**

United States (Headquarter)

6616 Owens Drive, Pleasanton, CA 94588, USA

Toll-Free (USA)  : +1 (844) 856-3350, customersupport@dbiosys.com



Contact your distributor
for more information
and visit us at www.dbiosys.com

Europe

Germany
Lake Constance

Cell: +49 7557 929 3915
wolfgang.vogel@dbiosys.com

Mexico

Avenida Aztecas # 95, Casa 4
Colonia Pueblo de los Reyes
Alcaldia Coyoacan
CDMX, C.P. 04330

Cell: +525535709693
a.herrera@dbiosys.com

India

612, Eden Square
St. John's Road
Secunderabad – 500003
Tel: 040-40078333

Cell: +91 9958293222
anandam@dbiosys.com

Canada

32 Lemsford Drive
Markham, ON L3S 4H4

Cell: +1 416 219 2035
hosna.mujadidi@dbiosys.com